



PRODUCT CATALOGUE

Trusted rehabilitation solutions
since 1976. Manufacturer of Class I
medical devices produced under
an ISO 13485 certified quality
management system.

WHO WE ARE

For 50 years, RMS (Rehabilitation Manufacturing Services) has built a reputation for dependable, clinically proven postural and mobility products, trusted across established healthcare organisations such as the NHS. As a British manufacturer, we maintain full control over design and production, ensuring consistent quality, reliable supply and the confidence that comes with proven performance.

Our focused range, including the Marley buggy, Kaye walker, Virtus and Caddy seating systems, alongside a comprehensive selection of postural support accessories, is designed to integrate seamlessly into existing portfolios.

Whether you are extending your offering or introducing proven products into a new market, RMS provides a flexible, responsive approach, supported by technical expertise, training and the ability to accommodate specific requirements where needed.

WHY RMS

- 50 years of proven manufacturing experience
- British-made, quality-controlled production
- Flexible approach to market and customer requirements
- ISO 13485 certified for quality and safety
- Training and technical support for partners


A STRONG PARTNERSHIP

RMS is more than a product supplier. With full UK based manufacturing and 50 years of experience, we provide a dependable and flexible foundation for businesses looking to strengthen their offering with proven postural and mobility solutions.

Our approach is built around consistency, responsiveness and collaboration, supporting partners with the confidence that comes from working with an established and trusted manufacturer.



ISO 13485 Certified



We operate under **ISO 13485** standards, ensuring every device meets strict medical and manufacturing requirements

- Extend your range of mobility and postural solutions without the need for product development
- Consistent, quality-controlled manufacturing in the UK
- Reliable supply chain from an experienced and established manufacturer
- Flexible approach to meet specific market and customer requirements



The Kaye Walker is **recommended by therapists**, trusted by the NHS and used in schools, clinics, homes and therapy centres around the world.

It fosters independence and participation in meaningful activities and supports users from early years to adulthood.




Trusted by clinicians for over 40 years

The Kaye Walker has been a **proven and trusted design** for over 40 years and it is now manufactured in the UK by RMS. It combines medical-grade production, real-world results and a legacy of clinical trust.

By tailoring the frame and support structures to each user's anthropometry, gait pattern and clinical needs, the Kaye Walker delivers a **uniquely effective geometry**.

Key Benefits

- **Supports upright posture:** the rear positioned frame encourages extension and natural trunk alignment
- **Built for real environments:** works in schools, clinics, therapy settings and home use
- **Customisable setup:** wheel kits, hand grips, forearm supports and hip guides to match user needs



The Kaye Walker is available in seven size options, each paired with a specific colour for easy identification:

W0

White

W1

Mint

W2

Orange

W3

Blue

W4L

Grey

W4H

Grey

W5

Aluminium



Adjustability & Growth


- Handle height adjustment
- Stability maintained at taller settings
- Tailoring frame and support structures to match the user's individual anthropometry, gait pattern and clinical needs

Frame & Structure

- Rear-positioned for upright gait
- Rigid aerospace-grade aluminium frame
- Seven sizes from paediatric to adult

Everyday Practicality

- Folds flat easily for transport/storage
- Sturdy, yet lightweight (from 3.0 kg)
- Easy-clean finish and low maintenance



	Height with Standard Wheels	Height with Activity Wheels	Hand Grip Width	FRAME		
				Width	Length	Weight
W0	37 - 46 cm	N/A	34 cm	58 - 60 cm	52 - 59 cm	3.00 kg
W1	41 - 55 cm	45 - 59 cm	34 cm	60 - 62 cm	56 - 62 cm	3.00 kg
W2	48 - 64 cm	52 - 68 cm	34 cm	58 - 60 cm	69 - 83 cm	3.85 kg
W3	59 - 78 cm	63 - 82 cm	38 cm	62 - 64 cm	77 - 89 cm	4.85 kg
W4L	72 - 92 cm	76 - 96 cm	40 cm	65 - 67 cm	84 - 97 cm	5.50 kg
W4H	81 - 100 cm	N/A	40 cm	67 - 69 cm	90 - 102 cm	7.00 kg
W5	90 - 105 cm	94 - 109 cm	47 cm	65 - 67 cm	102 - 108 cm	9.75 kg

USER	W0	W1	W2	W3	W4L	W4H	W5
Max Weight (kg)	27	27	39	60	82	98	113
Height (cm)	up to 95	91-122	107-137	129-152	150-179	150-185	175-192

► WHEEL OPTIONS

FRONT WHEELS



STANDARD 4" WHEELS

Direction Locking Swivel Castors



STANDARD 5" WHEELS

Forward and Back, Fixed



ACTIVITY 6" WHEELS

Direction Locking Swivel Castors



ACTIVITY 8" WHEELS

Forward and Back, Fixed

REAR WHEELS



STANDARD 5" WHEELS

Forward Only, Fixed



STANDARD 5" WHEELS

Forward and Back, Fixed



LEG FERRULE

Rear Conversion



ACTIVITY 6" WHEELS

Forward Only, Fixed

► ACCESSORIES



GUTTER ARMS

Provide upper-limb support for users who need extra help with balance or trunk stability while walking



HAND GRIPS

Provide a comfortable grip while promoting correct hand placement





HIP GUIDES

Support the hips to encourage correct alignment and improve sitting or standing posture



SEAT KIT & CUSHION

Secure rest platform
Folds away when walking



UNIVERSAL ANTI-TIP BARS

Prevent the frame from tipping backwards



REAR WHEEL DAMPENERS

Rear-wheel resistance to promote steady gait speed and improve stability while walking

▶ ACCESSORIES

ANTI PERCHING LOOP

Prevents the user from sitting on the frame, encouraging standing



SOFT SLING SUPPORT

Gives gentle support under the pelvis to help users move more smoothly across different floor surfaces

REMOVABLE CONTROL HANDLE

Allows a carer or therapist to guide and steady the walker, providing additional safety and reassurance during training or early walking stages





HOW TO ORDER

Kaye Walkers are available exclusively via prescription form, ensuring each walker is configured to meet the individual needs of the user. Please get in touch to request a copy.

[Download the Prescription Form](#)



marley

Rehabilitation Buggy



Built for Everyday Life

The Marley is an eye catching, practical, simple to use buggy that is designed to fit effortlessly into everyday family life.

With its fixed 90° hip angle, solid seat and backrest structure it ensures 90-90 seating is always maintained. Tilt in Space allows the user to be positioned anywhere within its 15° to 50° range without compromising their sitting position, whilst also aiding respiration, digestion and pressure relief.

The padded fabric covers are easily removable for washing and have been made from breathable spacer fabric that ensures the user always remains comfortable.

Key Features

- Adjustable seat pan
- Fixed 90° backrest angle
- 35° tilt in space
- Machine washable upholstery
- Padded pelvic belt
- Backrest extender
- Hip supports
- Upper trunk supports
- Shopping basket
- Height adjustable push handle

Transport & Storage

- All in one folding frame
- Quick release rear wheels

Options

- Available in two sizes
- Range of additional support accessories
- Rain and sun canopies

Compliance

- Fully tested to ISO7176
- RMS is ISO 13485 certified

► TECHNICAL SPECIFICATIONS

	SIZE 1		SIZE 2	
	Min	Max	Min	Max
Seat Width (with pelvic supports)	200 mm	270 mm	260 mm	320 mm
Seat Depth	230 mm	280 mm	290 mm	350 mm
Seat to Floor Height (upright)	490 mm		520 mm	
Backrest Height (from seat base)	510 mm		570 mm	
Backrest Height (with extender)	625 mm		720 mm	
Thoracic Width	120 mm	210 mm	170 mm	240 mm
Thoracic Height (from seat base)	150 mm	340 mm	160 mm	430 mm
Footplate Height	130 mm	280 mm	135 mm	330 mm
Footplate Depth	160 mm		180 mm	
Push Handle Height (from floor)	910 mm	1250 mm	925 mm	1270 mm
Tilt in Space Angle	15°	50°	15°	50°
Overall Width	580 mm		640 mm	
Overall Depth	790 mm		820 mm	
Overall Height	1070 mm		1160 mm	
Buggy Weight	15 kg		17 kg	
Maximum User Weight	35 kg		35 kg	
Buggy Folded Size	800 x 590 x 500 mm		830 x 640 x 500 mm	



Highly adjustable
for growth
and positioning




▶ ACCESSORIES

The Marley Buggy comes supplied as standard with:


- + Padded Pelvic Belt
- + Backrest Extender
- + Lateral Trunk Supports
- + Shopping Basket



Four Point Pelvic Belt
Side Pull | Metal Buckle



Dynamic Harness
Quick Release 20mm




Comfort Headrest
Height Adjustable



Leg Abductor



Back Insert



**Lateral Trunk
Supports**



Clear Tray



Grab Bar



Anklesures



Vent Tray Platform
Cargo Space



Vent Tray Platform
+ Oxygen Bottle Holder



Footplate Padding
25mm Foam

▶ ACCESSORIES



Sun Canopy



Rain Cover
(requires sun canopy)



Canopy + Rain Cover
Bundle

Reliable Postural Support for Everyday Mobility

POSTURAL CONTROL

As standard, the Marley is easily adjusted for correct postural support and growth. The seat is both width and depth adjustable incorporating hip guides whilst the backrest has fully adjustable trunk supports.

TILT IN SPACE

The Marley has a continuous tilt in space range of 15° to 50° to meet the individual requirements and comfort of the user. The tilt mechanism can be activated simply and smoothly.

Vary the seating angle throughout the day to support posture and function and to allow rest.



FOLDING FRAME

The Marley frame folds quickly and easily in one simple motion to a compact size, ideal for both transportation and storage.

The rear wheels can then be removed for an even smaller storage footprint, and to minimise weight.



HOW TO ORDER

Marley Buggies are available exclusively via prescription form, ensuring each one is configured to meet the individual needs of the user. Please get in touch to request a copy.

[Download the Prescription Form](#)

MY LYTE



Postural Support in a Lightweight Package

The My-Lyte Discovery is a hybrid buggy that combines the popular and proven **Shuttle Discovery** seating system with the My-Lyte outdoor pushchair base, a lightweight and highly durable frame manufactured by RMS in the UK.

Designed to support the demands of everyday family life, the My-Lyte Discovery provides reliable **postural support** in a practical and easy-to-use package. Offering some of the smallest sizes available for early intervention, it presents an attractive alternative to traditional special needs buggies without compromising on clinical support or functionality.

Its **modular design** can be adapted with a wide range of parts and accessories, making it easy to accommodate changing needs over time. Easy-access zips allow covers to be removed quickly for machine washing, while an optional wipe-clean fabric provides a practical solution for everyday use.


Key Features

- Fully supportive seat
- Adjustable modular design
- 90° to 160° combined range, with 40° tilt-in-space and 30° backrest recline
- Small achievable sizes to enable early intervention
- Comprehensive range of accessories
- Reversible seat

► TECHNICAL SPECIFICATIONS

The My-Lyte Discovery features a reversible seat unit and a compact folding design for convenient transport and storage.

It has also been crash tested for use during school transportation, providing added reassurance for everyday travel.



SEAT UNIT

	Min	Max
Seat Width	130 mm	270 mm
Seat Depth	130 mm	300 mm
Backrest Height	240 mm	390 mm
Seat to Head Height	450 mm	580 mm
Thoracic Width	120 mm	250 mm
Footplate Height	140 mm	280 mm

MY-LYTE BASE

	Min	Max
Client Weight	-	25 kg
Push Handle from Floor	900 mm	1270 mm
Tilt in Space with Seat	up to 40°	
Width x Depth	640 x 840 mm	
Frame Weight	12 kg	
Folded Size L x W x H	850 x 630 x 500 mm	





Extensive Adjustment for Individual Support



▶ ACCESSORIES




Standard Headrest




Occipital Headrest



Lateral Support



Pommel Abductor



Hip Guides




Adduction Thigh Guides



Tray Table



Padded Tray Top



Grab Bar



2 Point Belts



4 Point Belts



Harnesses


▶ ACCESSORIES



Foot Straps



Anklesures



Footplate Cover



Footplate Padding



Vent Tray Platform

Cargo Space



Vent Tray Platform

+ Oxygen Bottle Holder



Sun Canopy



Rain Cover

Requires Sun Canopy

The lightweight My-Lyte frame folds smoothly and fits into the boot of most cars, making everyday transport easier for parents and carers.

Download the Prescription Form



Buddy

CLASSROOM CHAIR



Designed for Learning, Built to Last

The Buddy Classroom Chair is an affordable seating solution designed for users requiring mild levels of postural support. Suitable for use in schools, therapy environments and at home, it combines practical adjustability, a wide range of accessories and durable construction in a cost-effective package.

Available in six sizes and compatible with furniture manufactured to BS EN 1729 standards, the Buddy integrates seamlessly into a variety of educational and care settings.

Key Features

- **Durable** polypropylene construction
- **Easy to clean** and maintain
- **Six sizes:** from early years to adulthood
- **Interchangeable** seat cushions
- Wide range of **accessories**
- **Affordable** seating solution
- **Highly recyclable**



Designed for everyday life

Manufactured from durable polypropylene, the Buddy is easy to clean, simple to maintain and built to withstand the demands of everyday use.

The Buddy has been developed to integrate seamlessly into classrooms, therapy spaces and home environments.

Thoughtful Colour Choices

Available in four carefully considered colours that balance modern classroom aesthetics with colours children naturally respond to. By moving away from a clinical appearance, the Buddy helps create a seating solution that children are happy to use as part of their everyday routine.



Mobility to suit your needs

The Buddy is available with either a static or mobile base, with the mobile version featuring two fixed silent front wheels and two rear swivel castors with lockable brakes for smooth manoeuvrability and easy positioning.

STATIC BASE



WHEELED BASE



	Seat Height	Seat Depth	Seat Width	Backrest Width	Backrest Height	Max User Weight
Size 1	21.5 - 26.5 cm	19 - 24.5 cm	24.5 cm	22 cm	20.5 cm	25 kg
Size 2	26 - 31.5 cm	25.5 - 31.5 cm	30 cm	27.5 cm	29.5 cm	45 kg
Size 3	31 - 39 cm	30 - 36.5 cm	34.5 cm	32 cm	35 cm	50 kg
Size 4	33.5 - 44 cm	34 - 41.5 cm	37.5 cm	35 cm	40 cm	60 kg
Size 5	38.5 - 47.5 cm	40 - 48.5 cm	43 cm	40.5 cm	47 cm	70 kg
Size 6	43 - 53 cm	46 - 55 cm	47.5 cm	45 cm	54 cm	80 kg

Flexible accessory options

A wide range of accessories allows the Buddy to be adapted to suit individual requirements. Integrated fixing points throughout the chair enable the addition of accessories such as headrests, armrests, chest straps, wheels, brakes and push handles, providing additional support and functionality where required.



Flat Cushion



Ramped Cushion



Armrest Infill Pads



Black Tray



Push Handle



Flip-Up Footrest



Stability Skis



Flip Down Pommel



Headrests



postureguard®
Complete range of belts
and harnesses



HOW TO ORDER

Buddy Chairs are available exclusively via prescription form, ensuring each one is configured to meet the individual needs of the user. Please get in touch to request a copy.

[Download the Prescription Form](#)



Made-to-Measure Support with Exceptional Adjustability

The Virtus Seating System is a highly adjustable, made-to-measure seating solution designed to support complex postural needs from early childhood through to adulthood. Its robust one-piece construction provides a stable foundation for individual positioning, with extensive adjustment to accommodate growth and changing requirements. Virtus is particularly suited to users who require high levels of postural support throughout daily activities.

A choice of fixed or dynamic backrest allows the system to be configured around the user's individual postural and movement needs. Its integrated interface system enables the seating unit to be transferred easily between compatible mobility bases, including powered wheelchairs, tilt-in-space bases and High-Low Bases, helping to maintain consistent support across different activities and environments.

Key Features

- Individually made to measure for each user
- Four sizes supporting users from around 1 year old through to adulthood
- Highly adjustable design to accommodate growth and changing postural needs
- Dynamic seating technology to support movement and tone management

Flexible Integration Across Multiple Mobility Bases

The Virtus Seating System can be paired with either the powered or manual High-Low Base, providing flexible positioning and support across a wide range of clinical and everyday environments. For additional versatility, Virtus can also be supplied with an interface that enables compatibility with virtually any suitable mobility base.

This adaptable approach allows the seating system to be transferred between multiple setups, including powered wheelchairs, tilt-in-space mobility bases and indoor or outdoor chassis systems, helping to maintain consistent postural support across different activities and environments.

Consistent
seating support
across multiple
environments

Max. Client Weight:

HI/LO BASE 90kg

MOBILITY BASE 135kg



Dynamic Support

The Virtus dynamic backrest is designed for users with high tone, involuntary movements or extensor patterns who may benefit from controlled movement whilst seated.

Its dynamic bearing system allows the backrest to respond to the user's movements, helping to manage the forces generated whilst maintaining postural support. This can improve comfort, positioning and overall seating tolerance.



Long-Term Adjustability

To support long-term use and future growth, the overall backrest height is individually configured using the user's shoulder height measurement with an additional 7.5cm growth allowance incorporated into the design.

	Size 1	Size 2	Size 3	Size 4
WHEELCHAIR SIZE	36	40	45	50
SEAT DEPTH	19 - 37 cm	26 - 45 cm	34 - 52 cm	34 - 52 cm
PELVIC WIDTH	21 - 32 cm	25 - 36 cm	30 - 41 cm	35 - 46 cm
THORACIC WIDTH	19 - 36 cm	25 - 40 cm	31 - 46 cm	35 - 51 cm
SACRAL PAD HEIGHT	10 cm	12.5 cm	15 cm	15 cm

Modular **Postural Support**

The Virtus Seating System is supplied with a comprehensive range of postural support accessories, enabling highly individual positioning and clinical configuration.

To accommodate the wide range of available options and customisations, Virtus seating systems are specified using our intuitive and interactive prescription form.

Trunk & Thoracic Support


- Fixed and swing-away laterals
- Shoulder and tricep supports
- Harness systems


Pelvic Positioning

- 2-point and 4-point belts
- Stabiliser belts
- Pelvic pads

Cushions & Covers

- Flat and contoured cushions
- Multiple foam options
- Washable and wipe-clean cover options





Head Support

- Integrated Headrest
- Modular headrests
- Chailey Rollercoaster
- Adjustable mounting

Arm Support & Trays

- Adjustable arm supports
- Padded and moulded tray options
- Functional support for everyday activities

Leg & Foot Support

- Pommels
- Knee abductors
- Dynamic footrests
- Anklesures and foot support options

Introduce flexible, clinically proven seating into your portfolio

Flexible Integration

- Powered wheelchairs
- Manual High-Low Bases
- Powered High-Low Bases
- Indoor chassis
- Outdoor mobility bases

Manufactured in the UK

- ISO 13485 certified manufacturing
- Quality-controlled production
- Durable components designed for long-term use
- Responsive technical support and aftersales assistance





HOW TO ORDER

Virtus Seating Systems are available exclusively via prescription form, ensuring each one is configured to meet the individual needs of the user. Please get in touch to request a copy.

[Download the Prescription Form](#)



VRTUS

⤴ HIGH - LOW ⤵



Powered Positioning with Everyday Practicality

The Virtus Powered High-Low Base is a robust and highly adaptable indoor mobility solution designed to support the demands of everyday home, school and clinical environments. Combining powered height adjustment and tilt functionality with a durable heavy-duty frame, the system provides flexible positioning whilst maintaining ease of use for carers and clinicians.

Designed as part of the wider Virtus seating range, the base integrates directly with the Virtus Seating System and can also be configured to support a wide range of third-party seating solutions. Its smooth manoeuvrability, substantial weight capacity and powered adjustment features help deliver practical, long-term support across a variety of care settings.

Key Features

- Powered height adjustment and tilt functionality
- Heavy-duty frame for everyday clinical and home use
- Direct Virtus Seating System integration
- Compatible with third-party seating systems
- Smooth manoeuvrability with lockable castors

▶ TECHNICAL SPECIFICATIONS

↖ ↗ **Base Size**
Width: 65 cm
↙ ↘ Depth: 83 cm



Weight Limit

User: Max. 90 kg
+ Seat: Max. 120 kg

Tilt-in-space

5° of Forward Tilt
25° of Backward Tilt

Height Adjustment

37.5cm at lowest
62.5cm at highest

25cm

20cm

Adjustable Push Handle

Min. 70cm
Max. 90cm

Castor Wheels

Lockable Ø125mm

Battery



- Up to 50 lift cycles per charge
- Can be operated whilst plugged in
- Integrated cable storage basket

POWERED ADJUSTMENT

The Virtus Powered High-Low Base provides smooth powered height adjustment with a 25cm lift range, supporting transfers, interaction and everyday positioning across home, school and clinical environments.

ENGINEERED FOR DAILY USE

The heavy-duty steel frame combines robust construction with smooth manoeuvrability, lockable castors and threshold clearance for practical indoor mobility.

FLEXIBLE SEATING INTEGRATION

The base interfaces directly with the Virtus Seating System and can also support a wide range of third-party seating solutions, providing a highly adaptable platform for different clinical and mobility requirements.

POWERED TILT

Powered forward and rear tilt provides up to 5° forward tilt and 25° rearward tilt, supporting positioning, comfort and hoisted transfers whilst helping maintain user stability.



Manual Base Available

User + Seat Max. 60kg

-10° - 30° Tilt

HOW TO ORDER

Virtus High-Low bases are available exclusively via prescription form, ensuring each one is configured to meet the individual needs of the user. Please get in touch to request a copy.

Download the Prescription Form



CADDY



Flexible Seating with Early Intervention Capability

The Caddy Seating System is a lightweight, highly modular seating solution designed to provide flexible postural support for a wide range of clinical needs. Its two-piece construction allows the backrest to be separated for easier transport, handling and storage, making it a practical solution for everyday use across home, school and community environments.

With sizes starting from as little as 20cm x 20cm, the Caddy is particularly well suited to early intervention and growing users. The system combines extensive adjustment, configurable support options and growth capability within a lighter and more transportable package, offering a flexible seating solution for users who do not require the advanced dynamic support features of the Virtus system.

Key Features

- Highly modular seating system with extensive adjustment
- Sizes from 20cm x 20cm for early intervention applications
- Growth adjustment built into seat depth and back height
- Wide range of configurable postural support options

Lightweight Two-Piece Design

The Caddy Seating System features a practical two-piece construction that allows the backrest to be separated from the seat base for easier transport, storage and handling.

Its lighter-weight design makes the system well suited to everyday use across home, school and community environments whilst maintaining stable and adaptable postural support.



Lightweight,
Modular and
Adaptable



Designed for Early Intervention

With sizes starting from as little as 20cm x 20cm, the Caddy is particularly suited to early intervention applications and smaller users. Its highly adjustable design enables clinicians and providers to achieve accurate positioning and support from an early stage whilst accommodating changing needs over time.

Cost-Effective Modular Seating

The Caddy Seating System provides a practical and cost-effective approach to postural seating without compromising on configurability or clinical support. Its modular design, lightweight construction and flexible specification make it an accessible solution for providers requiring adaptable seating across a wide range of user needs and environments.



Modular **Postural Support**

The Caddy Seating System supports a wide range of postural and positioning requirements through its extensive selection of configurable support options and accessories.


Typically fitted directly to manual wheelchair seat rails, the system can also be supplied with an interface and rail solution for selected powered wheelchair applications, providing additional flexibility across different mobility setups.


Pelvic Positioning

- Pelvic pads and supports
- Multiple bracket reductions
- Pelvic belts and stabilisers
- Configurable positioning options

Cushions & Covers

- Flat and contoured cushions
- Multiple foam options
- Washable and wipe-clean cover options





Head Support

- Integrated Headrest
- Modular headrests
- Chailey Rollercoaster
- Adjustable mounting

Trunk & Thoracic Support

- Fixed and swing-away laterals
- Flat, curved and wrap-around pads
- Shoulder and tricep supports
- Harness and upper body support

Leg & Lower Limb Support

- Leg abductors and pommels
- Knee adductors
- Anklesures

HOW TO ORDER


Caddy Seating Systems are available exclusively via prescription form, ensuring each one is configured to meet the individual needs of the user. Please get in touch to request a copy.

Download the Prescription Form

POSTURAL SUPPORT SOLUTIONS & ACCESSORIES


Developed through decades of in-house design and close collaboration with clinicians, the RMS range of postural support accessories provides flexible solutions for users with a wide variety of positioning needs. From simple positioning aids to highly specialised support systems, every component is designed to work together to achieve the most appropriate clinical outcome for each individual.

Our comprehensive range supports the body from head to feet, including headrests, trunk supports, pelvic positioning, harnesses, belts and lower limb supports. Compatible with a wide variety of seating systems and wheelchairs, RMS accessories can be specified individually or combined to create tailored postural support solutions.



Manufactured in the UK using quality controlled processes, they offer the reliability, flexibility and long-term performance trusted by clinicians and providers worldwide.

From Minimal to Advanced Support



▶ HEAD SUPPORT

Type A Headrest

Curved support providing medium to good head position. Can be fitted with one of our headrest covers and all types of headrest mounts.



Headrest Cover

SMALL L8-P18-S-CV

LARGE L8-P18-L-CV

SMALL L8-P18-S-#


LARGE L8-P18-L-#

A	B	C	D	E
210 mm	95 mm	95 mm	60 mm	25 mm
275 mm	120 mm	115 mm	70 mm	25 mm

Replace the # with the letter corresponding to the mount type on page 62.

Type D Headrest

Shaped headrest providing good comfort combined with minimal positioning. Can be fitted with comfortable headrest covers and all types of headrest mounts.



Headrest Cover

SMALL L8-P21-S-CV

LARGE L8-P21-L-CV


SMALL L8-P21-S-#

LARGE L8-P21-L-#

A	B	C	D	E
195 mm	50 mm	100 mm	60 mm	30 mm
300 mm	75 mm	150 mm	90 mm	45 mm

Replace the # with the letter corresponding to the mount type on page 62.

Type E Headrest



Contoured combination head and neck rest giving high levels of comfort and support.


All of our headrests are made from premium and robust materials - making them comfortable and highly functional.

ONE SIZE L8-P22-#

A	B	C	D	E
135 mm	75 mm	135 mm	70 mm	35 mm

Replace the # with the letter corresponding to the mount type on **page 62**.

Type K Headrest



Shaped headrest providing good comfort combined with minimal positioning.

Can be fitted with comfortable headrest covers and all types of headrest mounts.

SMALL L8-P29-S-#
LARGE L8-P29-L-#

A	B	C	D	E
140-240 mm	100-115 mm	100 mm	60 mm	30 mm
150-280 mm	100-130 mm	150 mm	90 mm	45 mm

Replace the # with the letter corresponding to the mount type on **page 62**.

Type J Headrest



Slightly curved oblong headrest offering a basic level of support.

135 x 75 x 35mm L8-P17-#

Replace the # with the letter corresponding to the mount type on **page 62**.

▶ HEAD SUPPORT

 ORDER BY
PRESCRIPTION

Type G Headrest

Also known as the Contour Headrest, it features multiple adjustable links that can be shaped to meet individual positioning requirements.

Available in a complete range of sizes, it can be configured symmetrically or offset to either side, with a choice of standard height or a double-height configuration for increased support.

Due to its high level of adjustability, we recommend using the **prescription form** when specifying a Type G headrest.

Popular Choices:

SMALL L8-P24-S-#

LARGE L8-P24-L-#

Replace the # with the letter corresponding to the mount type you need on **page 62**.



Replacement Pad

SMALL L8-P33-S

LARGE L8-P33-L



Type G Headband

SMALL L8-P30-S

LARGE L8-P30-L

Chailey Rollercoaster Headrest

The Chailey Rollercoaster was designed by Chailey Heritage Clinical Services, an NHS centre of excellence in the clinical management of children and young people with complex physical disabilities.

It provides a unique, highly effective and clinically proven head support, with multiple stem mounting options and easy hand wheel adjustments. The headrest can remain in place during transfers by raising the arms, supporting both ease of use and consistent positioning.

The system can be mounted to a wide range of seating systems using either 13 mm or 15 mm stems. It has been successfully fitted to systems including the Virtus, CAPS, R82 Xpanda and Leckey Advanced seat units, as well as Rea Assist and Chunc wheelchairs and Chailey toilet seats from Smirthwaite.



Chin Prompt

 ORDER BY
PRESCRIPTION

▶ HEAD SUPPORT

Standard Mounting Options

The headrests can be fitted to a wide range of seating systems using our versatile mounting options.

Replace the # in the **headrest product code** with the letter corresponding to the required mount type shown below.

Standard Mount



Angle Adjustable



Multi-Axis Cranked



Multi-Axis Straight



Custom Mounting Options

These mounts are suitable for particular wheelchairs or situations where the standard options don't fit, or a custom solution is required.

Rea Azalea Head Mount

13mm stem receiver slots into Rea Azalea wheelchair backrest

L4-P31-F3



Ball Mount with Collar

13mm square chrome bracketry, supplied with stainless steel mounting collar



200mm L8-P01-F3-A

300mm L8-P01-F3-B

90° Stem

13mm square chrome bracket with tri-knob



300mm L8-P03-F3-A

400mm L8-P03-F3-B

500mm L8-P03-F3-C

Swan Neck Stem

13mm square chrome bracket with tri-knob



300mm L8-P04-F3-A

400mm L8-P04-F3-B

Angle Adjustable

Supplied with swivel collar and 13mm or 16mm square bracketry



13mm G3-08-03

16mm G3-08-03-16

Extension Kit

A simple kit to extend angle adjustable stem by one or more 60mm pivot links



60mm G3-08-03-EXT

Square Height Adjust Collar



13mm L8-P08-F3

16mm L8-P08-F3B

Multi-Axis Cranked

Height adjustment collar, four point fixing and mounting plate



340mm L8-P36-F3

Multi-Axis Straight

Height adjustment collar, four point fixing and mounting plate



300mm L8-P37-F3

Tube Mount Clamp

25mm | 1"

Switchable 13 mm or 16 mm stem receiver clamp for standard 25 mm (1") wheelchair cross bars



L4-P27-F4

▶ HEADREST MOUNTING BRACKETS



2 Point Fixing

13mm G3-04-06

16mm G3-03-19



4 Point Fixing

13mm L3-P05-F3-A

16mm L3-P05-F3-B



Multi Point

13mm L4-P26-F5-A

16mm L4-P26-F5-B

Clamp on Fixing

Supplied with tri-knob and round clamp bolt - fitting finished in chrome.

CLAMP RECEIVER

19mm | 13mm L8-P27-F3-19/13

19mm | 16mm L8-P27-F3-19/16

22mm | 13mm L8-P27-F3-22/13

22mm | 16mm L8-P27-F3-22/16

25mm | 13mm L8-P27-F3-25/13

25mm | 16mm L8-P27-F3-25/16



Specialist Mounts

- Mounting plate allows fitting to a backrest shell or channel back system (e.g. Vtrak)
- Suitable for both 13mm and 16mm headrest stems



Headrest Mount

L4-P44-SSHM



Head + Harness Mount

L4-P44-SSSHM



Harness Mounting Frame

Headrest stem clamp for harness positioning on canvas back chairs


13mm Stem L4-P35-A

20mm Stem L4-P35-B

Quick Release Mounting Frame

Allows the attachment of our headrests onto folding wheelchairs with canvas backrests.

This headrest mounting bracket is designed to fit headrests with either 13mm or 16mm mounting brackets, plus it has mounting options for harnesses to be fixed in the correct position.



Full Kit

10" to 14" QRMF 1/10 KIT

15" to 20" QRMF 1/15 KIT

This quick release mounting kit has successfully passed a front impact test and is suitable for clients whilst in transportation in motor vehicles.

▶ CHEST SUPPORT STRAPS

Padded chest straps provide adjustable, targeted support for users who do not require a full harness. **A tri-glide fitting kit is included as standard.** Other fitting kit options are available to purchase separately.



XS 28 x 8cm

Buckle and Pull Tight Fastening
TS-001-PG

Hook and Loop Fastening
TS-005-PG



S 38 x 10cm

Buckle and Pull Tight Fastening
TS-002-PG

Hook and Loop Fastening
TS-006-PG



M 47 x 13cm

Buckle and Pull Tight Fastening
TS-003-PG

Hook and Loop Fastening
TS-007-PG



L 57 x 16cm

Buckle and Pull Tight Fastening
TS-004-PG

Hook and Loop Fastening
TS-008-PG

▶ LATERAL SUPPORT



Clamp on Swing Away Laterals

Designed to provide lateral trunk support for users with scoliosis when using a canvas-style backrest.

Fitted as a pair or individually, our laterals attach directly onto 22mm or 25mm backrest frame tubing.

All of our high quality lateral supports include a set of soft neoprene covers.

	LINK	PAD	SIDE
SINGLE	L4-P32-	-	-
LINK SIZE:			
25 mm = 025			
50 mm = 050			
75 mm = 075			
PAD SIZE:			
Small = S			
Large = L			
SIDE:			
Left = LH			
Right = RH			

	LINK	PAD	
PAIR	L4-P32-	-	-P
LINK SIZE:			
25 mm = 025			
50 mm = 050			
75 mm = 075			
PAD SIZE:			
Small = S			
Large = L			

Ex: L4-P32-025-S-P



Extension Links

+ sleeve & fittings

	EXTENSION	REPLACEMENT
25mm	L4-P32-EXT-25	L4-P32-LINK-25
50mm	L4-P32-EXT-50	L4-P32-LINK-50
75mm	L4-P32-EXT-75	L4-P32-LINK-75




TRUNK SUPPORT

Made with reinforced neoprene material, these harnesses provide secure and comfortable support to the upper body, encouraging stability in an upright seated posture.

Tri-glide fitting kit included as standard. Other fitting kit options are available to purchase separately.

Sternum Harness

Quick Release, 25mm



S

M

L


H012-PG H004-PG H005-PG

	A	B	C	D	E
S	270	225	130	60	130
M	330	255	155	60	150
L	460	260	200	60	200

All units are displayed in mm

Slimline Harness

Quick Release, 25mm



S

M

L

H014-PG H015-PG H016-PG


	A	B	C	D	E
S	360	230	170	30	140
M	420	250	190	30	160
L	470	260	200	30	180

All units are displayed in mm

Sternum and **Slimline** harnesses are ideal for users who don't require a larger surface area to be supported and are likely to overheat.

Shoulder Harness

Quick Release, 25mm



S

M

L

H006-PG

H007-PG


H008-PG

	A	B	C
S	300	95	45
M	330	115	65
L	430	145	75

All units are displayed in mm

Dynamic Harness

Quick Release, 25mm



S

M

L

H009-PG

H010-PG

H011-PG

	A	B	C	D
S	110	180	160	45
M	145	230	200	50
L	175	260	230	60

All units are displayed in mm


A **Shoulder Harness** is ideal for individuals with a high PEG feed placement, as it helps ensure the harness does not interfere with the PEG feed.

It also gives a firmer orientation of force to the client's shoulder for kyphotic presentations.

▶ TRUNK SUPPORT

Buggy Dynamic

Quick Release, 20mm



S

M

H009-B-PG


H010-B-PG

	A	B	C	D
S	110	180	160	45
M	145	230	200	50
L	175	260	230	60

All units are displayed in mm

Dynamic Harness

Fixed Strap, 25mm



S

M

L

H009-FS-PG

H010-FS-PG

H011-FS-PG

	A	B	C	D
S	110	180	160	45
M	145	230	200	50
L	175	260	230	60


All units are displayed in mm

Dynamic harnesses provide support for controlled movement while maintaining postural support. The top fixed straps provide extra comfort.

The **Buggy Dynamic** is specifically designed for buggy and pushchair use with top straps having quick release, anti slip adjustable buckles.

3 Point Y Harness

Right, Quick Release, 25mm




S M L
H017-PG H018-PG H019-PG

	A	B	C
S	250	300	190
M	300	350	220
L	360	430	240

All units are displayed in mm

3 Point Y Harness

Left, Quick Release, 25mm



S M L
H020-PG H021-PG H022-PG

	A	B	C
S	250	300	190
M	300	350	220
L	360	430	240

All units are displayed in mm

The **3 point harnesses** are designed to provide targeted postural control with directional support. Available in either left or right orientation.

UNDER NO CIRCUMSTANCES should these harnesses be used as the sole method of securing the user in their wheelchair during transportation in a motor vehicle.

▶ CONTOUR TRAYS



Our Contour Trays are perfect for wheelchairs and seating systems where a flat, smooth and secure work or stabilising surface is required for the user.

Easy to clean and disinfect, great for infection control between multiple users.

Available as a 6mm clear polypropylene tray or a partially or fully padded tray.



	A	B	C	D
SMALL	56cm	51cm	33cm	15cm
MEDIUM	56cm	56cm	36cm	20cm
LARGE	61cm	61cm	43cm	22.5cm

Available with **three mounting options:**



DROP IN



SLIDE IN



SWING OVER

Bespoke solutions and **custom sizing**, manufactured in-house



▶ PELVIC SUPPORT

Maintaining neutral pelvic alignment is key to effective postural control. Pelvic positioning belts help reduce forward sliding and support functional seating outcomes.

Please use the guide below to select the appropriate belt, buckle and fitting kit.

01 Select your belt type



2-POINT PELVIC BELT



4-POINT PELVIC BELT



PELVIC STABILISER BELT

02 Select buckle type



STANDARD NON SLIP ANTI TAMPER

PLASTIC BUCKLE

Lightweight and rust-free, suited for outdoor or wet environments



METAL BUCKLE

Strong and durable, ideal for high-load applications

03 Select size

Please refer to the relevant detailed size charts in this brochure, on the prescription form or on our website.

04 Choose the fitting kit

- A STANDARD (Included)**
TRI-GLIDE: Simple and adjustable
- B FRAME CLAMPING KIT**
Attaches to the frame
- C TRIPLE SLOT MOUNTING PLATES**
Similar to the tri-glide but they attach to the clamping kit.
- D LOCKING CLAMPS**
Fast adjustment and secure lock, also attaches to the clamping kit.

▶ TWO POINT PELVIC BELT

A comfortable, padded pelvic belt providing secure support with straightforward clinical adjustment for users not needing a four point belt.

Pelvic Belt with Plastic Buckle

Padded | Dual Pull | Two Point

BPB-004-S-PG | BPB-004-M-PG | BPB-004-L-PG



Pelvic Belt with Metal Buckle

Padded | Side Pull | Two Point

BMB-001-S-PG | BMB-001-M-PG | BMB-001-L-PG



Pelvic Belt with Plastic Buckle

Padded | Side Pull | Two Point

BAT-007-S-PG | BAT-012-M-PG | BAT-010-L-PG



	Padding (Individual)	Webbing & Buckle
S	150 x 40mm	25mm
M	170 x 60mm	40mm
L	230 x 70mm	40mm

Pelvic Belt with Plastic Buckle

Padded | Rear Pull | Two Point

BPB-006-S-PG | BPB-006-M-PG | BPB-006-L-PG



Pelvic Belt with Plastic Buckle

Padded | Side Pull | Two Point

BPB-013-S-PG



Pelvic Belt with Metal Buckle

No Padding | Side Pull | Two Point

BMB-003-S-PG | BMB-003-M-PG | BMB-003-L-PG



► FOUR POINT PELVIC BELT

Provides four-point fixation to help maintain a neutral pelvic position and reduce forward sliding, with padded supports for improved comfort and pressure distribution.

Pelvic Belt with Plastic Buckle

Padded | Dual Pull | Four Point

[BPB-005-S-PG](#) | [BPB-005-M-PG](#) | [BPB-005-L-PG](#)



Pelvic Belt with Metal Buckle

Padded | Side Pull | Four Point

[BMB-002-S-PG](#) | [BMB-002-M-PG](#) | [BMB-002-L-PG](#)



	Padding (Individual)	Webbing & Buckle
S	150 x 40mm	25mm
M	170 x 60mm	40mm
L	230 x 70mm	40mm

All lower straps are 25mm wide

Pelvic Belt with Plastic Buckle

Padded | Side Pull | Four Point

BPB-014-S-PG



Pelvic Belt with Plastic Buckle

Padded | Side Pull | Four Point

BAT-008-S-PG | BAT-013-M-PG | BAT-011-L-PG




▶ PELVIC STABILISER BELT

Designed to deliver increased pelvic coverage and high-level support to stabilise the pelvis for more complex postural needs.

- A larger contact area helps stabilise the pelvis and provides greater control for users with stronger movement patterns.
- Provides stable, consistent support.
- Ideal for users with complex postural needs and users who struggle with effective positioning while using a two-point belt.

Pelvic Stabiliser Belt

20mm Buckle | Side Pull | Four Point




XXS

	Padding (Individual)	Webbing	Buckle
BPB-004-XXS-PG	125 x 98mm	20mm	20mm

Lower straps are 20mm wide


Pelvic Belt with Metal Buckle
25mm Buckle | Side Pull | Four Point



	Padding (Individual)	Webbing	Buckle
BMB-004-XS-PG	175 x 128mm	25mm	25mm
BMB-004-S-PG	195 x 158mm	25mm	25mm
BMB-004-M-25-40-PG	205 x 188mm	40mm	25mm
BMB-004-L-25-40-PG	225 x 220mm	40mm	25mm

Lower straps are 25mm wide

Pelvic Belt with Metal Buckle
40mm Buckle | Side Pull | Four Point



	Padding (Individual)	Webbing	Buckle
BMB-004-M-PG	205 x 188mm	40mm	40mm
BMB-004-L-PG	225 x 220mm	40mm	40mm

Lower straps are 25mm wide

▶ LIMB SUPPORT



Flip Down Pommel

The RMS Flip Down Pommel is designed to reduce leg adduction, help maintain a correct midline posture and prevent wind sweeping.

The pommel has a quick release button that flips it down and safely out of the way for quick transfers, whilst links allow for tooled adjustment.

S L3-P50-S **M** L3-P50-M **L** L3-P50-L



Slide in Pommel

Perfect for specialist seating, with a comfort block that can be moved in and out from the seat cushion, further helping to prevent leg sweeping and adduction.

	Length	Width	Height
S	9.3 cm	6.6 cm	9 cm
M	11.3 cm	8.2 cm	10 cm
L	12.8 cm	9.1 cm	10 cm

Slide in Short (8-10 cm)

S L3-P67-S-S **M** L3-P67-S-M **L** L3-P67-S-L

Slide in Tall (11-16 cm)

S L3-P67-L-S **M** L3-P67-L-M **L** L3-P67-L-L

Adductor Pads

Designed to effectively position the user's legs when seated in a wheelchair. The clamp-on mounting brackets fit a variety of tube sizes, and are fully adjustable, ensuring the pad is positioned at the perfect knee height for optimum support.



PAD TYPE	SIZE	OFFSET	CLAMP
EVAZOTE	L7-P01-		-

GEL	L7-P02-		-
------------	---------	----------------------	------------------------

PAD SIZE:	OFFSET:	CLAMP SIZE:
Small = S	25 mm = 25	19 mm = 19
Large = L	50 mm = 50	22 mm = 22
	75 mm = 75	25 mm = 25

Ex: L7-P01-S25-19

Stump Boards

Can be fitted to either left or right wheelchair footrest hangers and come with a black polyurethane flip-up pad. Supplied in a variety of angled stems from 0° to 30°, at 5° stepped increments with clamp sizes to suit almost all applications.

When ordering please state the hanger angle and stem size that you require so that our customer service team can provide you with the full part numbers and quotation.



L6-P03-#

► VIRTUS CUSHIONS

Complete seating solutions to accommodate a wide variety of postural requirements. From cost-effective flat cushions through to advanced adjustable contoured designs, the range provides flexibility to balance positioning needs, functionality and budget.

Options are available to help prevent forward sliding, provide additional pelvic and leg support and allow the cushion shape to be tailored to the individual user, delivering the appropriate combination of comfort, stability and postural support.



Flat Seat


L3-P52-#



Anti-Thrust

Adjustable

L3-P53-#



Contour

Adjustable

L3-P54-#


▶ DROP IN SEAT

Drop in seats are fitted as an alternative to a wheelchair seat canvas. They are mounted using four drop in brackets and secured using either a webbing strap or Lock & Latch interface brackets. They also allow the fitting of alternative seat cushions or accessories like a pommel.

Drop in Seat

With Hook and Loop Strips
and Drop Bracketry

L3-P16-#



▶ LOCK & LATCH

The Lock and Latch system enables Virtus and Caddy seating units to be quickly and securely attached to most standard wheelchair frames.

Designed to retain seat hooks in position, it features a simple push-pull locking pin mechanism that combines strength, durability and ease of use.



For 19mm Tubing LL19

For 22mm Tubing LL22

For 25mm Tubing LL25

▶ LIMB SUPPORT

We manufacture a range of upper and lower limb support solutions designed to promote optimum seating outcomes. These supports help maintain comfortable and appropriate limb positioning to enhance overall postural stability.

They are not intended for use as restraint systems.

Padded Footstraps (pair)

Two Piece Hook And Loop

Fitting Kit included



FS004S

S 135mm Padding

FS004L

L 160mm Padding

D'Ring Footstraps

Two Piece Hook And Loop

FS002 | 260 x 25mm



One Piece Hook and Loop

FS003S | 380 x 25mm

FS003L | 420 x 25mm



Anklesures (pair)



		Buckle Closure	Hook & Loop
XS	ø 49-64mm	AS001B	AS001
S	ø 51-68mm	AS002B	AS002
M	ø 56-75mm	AS003B	AS003
L	ø 75-94mm	AS004B	AS004
XL	ø 94-110mm	AS005B	AS005

Leg Abductor (pair)

Thigh Straps



	LAS-001
S	230mm Padding
	LAS-002
M	320mm Padding
	LAS-003
L	370mm Padding

Neoprene Calf Panel

One Piece Hook And Loop



	CS003
XS	75mm x 330mm
	CS004
S	100mm x 430mm
	CS005
M	150mm x 560mm
	CS006
L	150mm x 640mm

▶ FOOTPLATES


Designed to replace unsuitable stock wheelchair footplates, these multi-adjustable options provide the secure attachment of various foot support equipment. They are made from 12mm textured polypropylene to prevent slipping and feature pre-formed slots and holes for heel loops, foot straps and ankle guides.

All kits include adjustable mountings compatible with fixed or angle-adjustable stems. They can fit directly onto existing wheelchair stems, or we can supply new stems when required.

Slotted Footrests (Pair)

SMALL 225 x 115 mm MPF001/S-BPP

LARGE 275 x 130 mm MPF001/L-BPP



Adjustable Stems (Pair)

MPF002/16

Fixed Stems (Pair)

MPF003/16

▶ FOOTBOARDS

Designed to attach directly to the footplates to provide additional foot room, these footboards are manufactured from 12mm thick black polypropylene. They come complete with fixings and depth adjustment and are available as a single one-piece board or as a pair of individual ones.



One Piece Footboard

31 x 20 cm L6-P05-13/20-BPP

31 x 25 cm L6-P05-13/25-BPP

31 x 30 cm L6-P05-13/30-BPP

36 x 20 cm L6-P05-15/20-BPP

36 x 25 cm L6-P05-15/25-BPP

36 x 30 cm L6-P05-15/30-BPP

41 x 20 cm L6-P05-17/20-BPP

41 x 25 cm L6-P05-17/25-BPP

41 x 30 cm L6-P05-17/30-BPP



Individual Footboards

12.7 x 20 cm L6-P06-13/20-BPP

12.7 x 25 cm L6-P06-13/25-BPP

12.7 x 30 cm L6-P06-13/30-BPP

16.5 x 20 cm L6-P06-15/20-BPP

16.5 x 25 cm L6-P06-15/25-BPP

16.5 x 30 cm L6-P06-15/30-BPP

19 x 20 cm L6-P06-17/20-BPP

19 x 25 cm L6-P06-17/25-BPP

19 x 30 cm L6-P06-17/30-BPP



Dynamic Footboard

L6-P16-#



Comfort Beadbagg

L6-P14-#

► FOOTPLATE CUSTOMISATION



Footrest Raiser Blocks

16mm / 19mm L6-P02-F4-A

16mm / 22mm L6-P02-F4-B

16mm / 25mm L6-P02-F4-C

19mm / 19mm L6-P02-F4-D

19mm / 22mm L6-P02-F4-E


19mm / 25mm L6-P02-F4-F

22mm / 22mm L6-P02-F4-G

22mm / 25mm L6-P02-F4-H

16mm / 16mm L6-P02-F4-I

25mm / 25mm L6-P02-F4-J



Footrest Stems

L6-P01-F3

Supplied in pairs to fit 19mm hangers or our Raiser blocks to fit other sizes.



► CYLINDER CARRIERS

Our universal wheelchair cylinder carriers are compatible with BOC (DD, CD & PD), Medigas (PD & MB), Air Products (PA2) cylinders, plus all other medical oxygen cylinders up to 535mm in height.



**Cylinder Carrier
for Tilt-in-space
Wheelchairs**

QR-OX-C

Suitable for most tilt in space wheelchairs, incorporating additional bracing straps and a depth-adjustable front enclosure panel.

The carrier allows the bottle to be in a fully controlled position even with the wheelchair in a reclined position.



**Cylinder Carrier
for Standard
Wheelchairs**

QR-OX-A

Designed to hold a wide range of medical oxygen cylinders.

This quick release oxygen cylinder carrier will fit most standard non-tilt NHS wheelchairs, and is easily removable for quick folding and storage.

Your Brand. Our Manufacturing.

Many of our customers know RMS for manufacturing high-quality rehabilitation products. What is less well known is that we can also manufacture products under your brand.

Whether you are looking to strengthen your own product range, create a more consistent brand identity or offer your customers exclusive products, our branding and OEM manufacturing service provides a simple and cost-effective solution.



We already manufacture products for a number of well-known companies, allowing them to sell proven products under their own brand, supported by our manufacturing expertise and quality systems.

Our service can include:

- Your company branding and product labelling
- Custom packaging
- Product literature and instructions
- OEM manufacturing
- Bespoke product development
- Flexible production volumes, subject to minimum order quantities (MOQs)

From a single product line to a complete product range, we will work with you to understand your requirements and develop a solution that fits your business. Projects are subject to **minimum order quantities**, depending on the product and level of customisation required.

With over 50 years of manufacturing experience, ISO 13485 certification and customers in many countries around the world, you can be confident that your brand is supported by a trusted manufacturing partner.

Let's Discuss Your Requirements

Ready to Work Together?

From a single product to a complete branded range, RMS provides the manufacturing expertise, technical support and flexibility to help your business grow. We'd be delighted to discuss your requirements and explore how we can work together.

Get in touch with our team today.

 +44 (0) 1795 477 280

 www.rms-rehab.co.uk

 sales@rms-rehab.co.uk



TABLE OF CONTENTS

COMPLETE PRODUCTS

• Kaye Walker	6
• Marley Buggy	16
• My-Lyte Buggy	24
• Buddy Classroom Chair	32
• Virtus Seating System	38
• Virtus High-Low Base	46
• Caddy Seating System	50

POSTURAL SUPPORT ACCESSORIES

Head Support

• Headrests Type A, D, E, K, J	58
• Type G Headrest	60
• Rollercoaster Headrest	61
• Headrest Mounting Options	62

Trunk Support

• Chest Straps	66
• Lateral Support	67
• Harnesses	68

Trays	72
--------------	----

Pelvic Support	73
-----------------------	----

• Two Point Belts	74
• Four Point Belts	76
• Pelvic Stabiliser Belts	78

Limb Support

• Pommels	80
• Adductor Pads	81
• Stump Boards	81
• Virtus Cushions	82
• Drop In Seat	83
• Footstraps	84
• Anklesures	85
• Footplates	86
• Footboards	87
• Cylinder Carriers	89

YEARS OF
50
IMPROVING LIVES

RMS
est.1976

@RMSRehab



01795 477 280



www.rms-rehab.co.uk



sales@rms-rehab.co.uk



Thompson House, Unit 10 Styles Close,
Sittingbourne, Kent, ME10 3BF