

PRODUCT CATALOGUE

Postural support solutions engineered for
safety, comfort and long-term reliability.
Designed and delivered by experts.



posture
guard®


Clinical Benefits

The PostureGuard series is built on **decades of clinical experience** and practical collaboration. Over the past **50 years**, RMS's highly skilled team has continually refined and advanced the product range, drawing upon both internal expertise and valuable feedback from clinicians working within Wheelchair Services across the United Kingdom.

RMS Expertise and Collaborative Development

In recent years, the PostureGuard series has been further refined through **close collaboration** with experienced paediatric clinicians working in the field of posture and mobility. This has been achieved via dedicated engagement within schools and the broader community, ensuring the range remains responsive to the evolving and diverse needs of its users.

Key Clinical Considerations



PELVIC STABILITY

The pelvis forms the foundation for optimal postural control. The primary objective is to stabilise the pelvis and, wherever possible, maintain neutral pelvic alignment.

Selecting an appropriate pelvic positioning belt can help prevent or reduce posterior pelvic tilt, which in turn minimises sliding forwards and promotes a safer, more functional seated posture.

OPTIMUM TRUNK ALIGNMENT

Postural belts and harnesses can be prescribed to support an upright, aligned posture. These accessories are often employed alongside modular seating systems, supportive backrests and contoured seat cushions.

Supporting the trunk in its optimal alignment enhances the functional potential of the upper limbs and head.

SUPPORTING FUNCTION AND INDEPENDENCE

Optimal trunk alignment and a secure and stable pelvis can improve functional outcomes.

Providing proximal stability via appropriate support and alignment enhances distal stability, facilitating everyday activities.

CREATING STABILITY

The correct selection of postural belts and harnesses can reduce the risk of sliding forwards while seated, thereby helping to protect the skin from shear forces.

The best results can be achieved when the correct accessories are used in combination with appropriate seating solutions.

SUPPORTING RESPIRATORY FUNCTION

Thoracic (lateral) supports and the appropriate postural harnesses can help maximise lung capacity by relieving pressure on the internal respiratory organs, thereby supporting optimal respiratory function.



Did you know?

As a British manufacturer, we can produce bespoke belts and harnesses, or provide modifications, to meet individual requirements.

Get in touch to find out more.



“

'RMS is committed to ensuring that no shortcuts are taken when providing the best postural outcomes for customers.'

Our passion lies in creating solutions that deliver truly customised results, especially for those with individual and complex needs.'

Alan Wombell, Director

”



Scan the QR code to view the full RMS range online, including fitting guides, technical specifications and support resources.



Strength in **every stitch**



Choosing the right product helps optimise fit, ease of use and clinical effectiveness

In postural support, failure is not always obvious. It may be as subtle as a slipping belt, a rotating harness or a loosening joint.

That is why every PostureGuard product is engineered to maintain its shape, position and performance over time, drawing on **over 50 years of RMS manufacturing expertise** and a clear understanding of where stresses occur.

At our UK facility, each seam is manufactured for strength where it matters most.

Bar tacking and box stitching reinforce high-load points such as belt ends, harness anchors and buckle junctions.

RMS operates a Quality Management System certified to **ISO 13485:2016**, covering the design, manufacture, supply and aftercare of our rehabilitation and postural support equipment – giving clinicians and occupational therapists confidence in consistent, high-quality outcomes.

It's not about looking strong.
It's about **performing** under load, day after day.

- UK-made in our own workshop for complete quality control
- Reinforced at every point of clinical load
- Designed to maintain alignment and safety under tension
- Over 50 years of RMS manufacturing experience



PELVIC BELTS

A secure pelvic positioning belt is fundamental to functional seating outcomes. Without a securely supported pelvis optimum head, trunk and limb control is compromised.

All pelvic belts in this brochure are part of the PostureGuard range and are **manufactured in the UK** by RMS. They are designed with comfort, durability and ease of clinical adjustment as the key priorities.

Quick ordering decision guide

Follow the steps below to identify the most appropriate PostureGuard pelvic belt and corresponding product code.

01 Select your belt type



2-POINT PELVIC BELT



4-POINT PELVIC BELT



PELVIC STABILISER BELT

02 Select buckle type



STANDARD



NON SLIP



ANTI TAMPER

PLASTIC BUCKLE

Lightweight and rust-free, suited for outdoor or wet environments



METAL BUCKLE

Strong and durable, ideal for high-load applications

03 Select size

Please refer to the relevant detailed size charts in this brochure, on the prescription form or on our website.

04 Choose the fitting kit

- A STANDARD (Included)**
TRI-GLIDE: Simple and adjustable
- B FRAME CLAMPING KIT**
Attaches to the frame
- C TRIPLE SLOT MOUNTING PLATES**
Similar to the tri-glide but they attach to the clamping kit.
- D LOCKING CLAMPS**
Fast adjustment and secure lock, also attaches to the clamping kit.

▶ TWO POINT PELVIC BELT

A comfortable, padded pelvic belt providing secure support with straightforward clinical adjustment for users not needing a four point belt.

Two-point pelvic belts are available with a range of buckle and adjustment options, allowing selection to be matched to the user, setting and clinical requirements.

	Padding (Individual)	Webbing & Buckle
S	150 x 40mm	25mm
M	170 x 60mm	40mm
L	230 x 70mm	40mm

Pelvic Belt with Plastic Buckle

Padded | Dual Pull | Two Point

BPB-004-**S**-PG | BPB-004-**M**-PG | BPB-004-**L**-PG



Pelvic Belt with Metal Buckle

Padded | Side Pull | Two Point

BMB-001-**S**-PG | BMB-001-**M**-PG | BMB-001-**L**-PG



Pelvic Belt with Plastic Buckle

Padded | Side Pull | Two Point

BAT-007-**S**-PG | BAT-012-**M**-PG | BAT-010-**L**-PG



Pelvic Belt with Plastic Buckle

Padded | Side Pull | Two Point

BPB-013-**S**-PG



Pelvic Belt with Plastic Buckle

Padded | Rear Pull | Two Point

BPB-006-**S**-PG | BPB-006-**M**-PG | BPB-006-**L**-PG



Pelvic Belt with Metal Buckle

No Padding | Side Pull | Two Point

BMB-003-**S**-PG | BMB-003-**M**-PG | BMB-003-**L**-PG



Tri-Glide Fitting Kit included

▶ FOUR POINT PELVIC BELT

Provides four-point fixation to support the pelvis in a neutral, upright position to help prevent sliding forwards for more complex postural needs.

Four anchor points for enhanced pelvic control, with padding to enhance comfort and reduce risk of localised pressure on bony prominence.

	Padding (Individual)	Webbing & Buckle
S	150 x 40mm	25mm
M	170 x 60mm	40mm
L	230 x 70mm	40mm

All lower straps are 25mm wide

Pelvic Belt with Plastic Buckle

Padded | Dual Pull | Four Point

BPB-005-**S**-PG | BPB-005-**M**-PG | BPB-005-**L**-PG



Pelvic Belt with Metal Buckle

Padded | Side Pull | Four Point

BMB-002-**S**-PG | BMB-002-**M**-PG | BMB-002-**L**-PG



Pelvic Belt with Plastic Buckle

Padded | Side Pull | Four Point

BPB-014-**S**-PG



Pelvic Belt with Plastic Buckle

Padded | Side Pull | Four Point

BAT-008-**S**-PG | BAT-013-**M**-PG | BAT-011-**L**-PG



Tri-Glide Fitting Kit included

▶ PELVIC STABILISER BELT

Designed to deliver increased pelvic coverage and high-level support to stabilise the pelvis for more complex postural needs.

- Larger surface coverage pelvic belt designed with increased surface area to stabilise the pelvis and provide greater contact and control around the pelvis for those with stronger movement patterns.
- Four anchor points provide stable, consistent support.
- Ideal for users with complex postural needs and users who struggle with effective positioning while using a two-point belt.

Pelvic Stabiliser Belt

20mm Buckle | Side Pull | Four Point



XXS

BPB-004-XXS-PG

Padding
(Individual)

125 x 98mm

Webbing

20mm

Buckle

20mm

Lower straps are 20mm wide

Pelvic Stabiliser Belt

25mm Buckle | Side Pull | Four Point



XS

BMB-004-XS-PG

Padding
(Individual)

175 x 128mm

Webbing

25mm

Buckle

25mm

S

BMB-004-S-PG

195 x 158mm

25mm

25mm

M

BMB-004-M-25-40-PG

205 x 188mm

40mm

25mm

L

BMB-004-L-25-40-PG

225 x 220mm

40mm

25mm

Lower straps are 25mm wide

Pelvic Stabiliser Belt

40mm Buckle | Side Pull | Four Point



M

BMB-004-M-PG

Padding
(Individual)

205 x 188mm

Webbing

40mm

Buckle

40mm

L

BMB-004-L-PG

225 x 220mm

40mm

40mm

Lower straps are 25mm wide

Tri-Glide Fitting Kit included

► DYNAMIC GROIN HARNESS

Groin Straps

Padded | Plastic Buckles

GH-001-PG | GH-002-PG | GH-003-PG | GH-004-PG



S	GH-001-PG 300 x 250mm
M	GH-002-PG 400 x 300mm
L	GH-003-PG 600 x 400mm
XL	GH-004-PG 800 x 500mm

PostureGuard Groin Straps are supplied with **a full fitting kit** as standard, these include adjustable locking clamps, triple mounting plates and frame clamps.

Designed to support pelvic positioning and reduce forward sliding through controlled support.

A groin harness should be the last resort and should only be considered when all other postural accessories have been trialled and considered.

► CHEST SUPPORT STRAPS

Padded chest straps provide adjustable, targeted support for users who do not require a full harness. Available with buckle and pull-tight fastening or hook and loop fastening.

Support Strap | Buckle and Pull Tight Fastening



XS 28 x 8cm TS-001-PG



S 38 x 10cm TS-002-PG



M 47 x 13cm TS-003-PG



L 57 x 16cm TS-004-PG

Support Strap | Hook and Loop Fastening



XS 28 x 8cm TS-005-PG



S 38 x 10cm TS-006-PG



M 47 x 13cm TS-007-PG



L 57 x 16cm TS-008-PG

Tri-Glide Fitting Kit included


► UPPER TRUNK BODY SUPPORT

We manufacture a range of upper trunk support harnesses designed to promote optimum seating outcomes. These harnesses provide secure and comfortable support to the upper body, encouraging stability in an upright seated posture.

Made with reinforced neoprene material for comfort and breathability.

Sternum Harness

Quick Release, 25mm




H012-PG H004-PG H005-PG

	A	B	C	D	E
S	270	225	130	60	130
M	330	255	155	60	150
L	460	260	200	60	200

All units are displayed in mm

Slimline Harness

Quick Release, 25mm




H014-PG H015-PG H016-PG

	A	B	C	D	E
S	360	230	170	30	140
M	420	250	190	30	160
L	470	260	200	30	180

All units are displayed in mm

Shoulder Harness

Quick Release, 25mm




H006-PG H007-PG H008-PG

	A	B	C
S	300	95	45
M	330	115	65
L	430	145	75

All units are displayed in mm

Dynamic Harness

Quick Release, 25mm




H009-PG H010-PG H011-PG

	A	B	C	D
S	110	180	160	45
M	145	230	200	50
L	175	260	230	60

All units are displayed in mm

Dynamic Harness

Fixed Strap, 25mm




H009-FS-PG H010-FS-PG H011-FS-PG

	A	B	C	D
S	110	180	160	45
M	145	230	200	50
L	175	260	230	60

All units are displayed in mm

Buggy Dynamic

Quick Release, 20mm



H009-B-PG H010-B-PG


	A	B	C	D
S	110	180	160	45
M	145	230	200	50
L	175	260	230	60

All units are displayed in mm

Tri-Glide Fitting Kit included

3 Point Y Harness

Right, Quick Release, 25mm



S

M

L


H017-PG H018-PG H019-PG

	A	B	C
S	250	300	190
M	300	350	220
L	360	430	240

All units are displayed in mm

3 Point Y Harness

Left, Quick Release, 25mm



S

M

L


H020-PG H021-PG H022-PG

	A	B	C
S	250	300	190
M	300	350	220
L	360	430	240

All units are displayed in mm

Mounting Frames

For Headrest Stems



13mm

L4-P35-A

20mm

L4-P35-B

Handicare

L4-P35-C

Sunrise Medical

L4-P35-D

Fits on headrest stems and provides upper mounting points for harnesses when the seating system doesn't allow it.

RECOMMENDATIONS

- **Sternum** and **Slimline** harnesses are ideal for users who don't require a larger surface area to be supported and are **likely to overheat**.
- A **Shoulder Harness** is ideal for individuals with a **high PEG feed placement**, as it helps ensure the harness does not interfere with the PEG feed. It also gives a firmer orientation of force to client's shoulder for kyphotic presentations.
- **Dynamic** harnesses provide support for **controlled movement** while maintaining postural support. The top **fixed straps** provide extra comfort.
- The **Buggy Dynamic** is specifically designed for buggy and pushchair use with top straps having quick release, anti slip adjustable buckles.
- The **3 point harnesses** are designed to provide targeted postural control with directional support. Available in either left or right orientation.

Use these simple steps to identify the most appropriate upper trunk harness.

- Select harness **type**
- Select the appropriate **size** by using the charts available for each model
- Order a suitable **fitting kit** if additional fixing options are required beyond the supplied Tri Glide system.
- Order by referencing the **product codes**

IMPORTANT PLEASE READ

UNDER NO CIRCUMSTANCES should these harness be used as the sole method of securing the user in their wheelchair during transportation in a motor vehicle.

► LIMB SUPPORT

We manufacture a range of upper and lower limb support solutions designed to promote optimum seating outcomes. These supports help maintain comfortable and appropriate limb positioning to enhance overall postural stability.

They are not intended for use as restraint systems.

Padded Footstraps (pair)

Two Piece Hook And Loop



Fitting Kit included

Used to support foot positioning and alignment within wheelchair and seating systems.

- S** FS004S-PG 135mm Padding
- L** FS004L-PG 160mm Padding

Neoprene Calf Panel

One Piece Hook And Loop



Provides lower leg support and positioning.

- XS** CS003 75mm x 330mm
- S** CS004 100mm x 430mm
- M** CS005 150mm x 560mm
- L** CS006 150mm x 640mm

D'Ring Footstraps

Two Piece Hook and Loop




Fitting Kit included

- FS002** 260 x 25mm Webbing

D'Ring Footstraps

One Piece Hook and Loop



Fitting Kit included

- S** FS003S 380 x 25mm
- L** FS003L 420 x 25mm

Leg Abductor (pair)

Thigh Straps



Fitting Kit included

Used to position the thigh and assist with hip abduction in conjunction with an appropriate postural support system.

- S** LAS-001-PG 230mm Padding
- M** LAS-002-PG 320mm Padding
- L** LAS-003-PG 370mm Padding

Houdini Strap (pair)

Adjustable Arm Straps



- HOU1-PG** 220mm Padding
- Used to keep the users' arms from extending outside the wheelchair, especially useful in transport.

Anklesures (pair) + Fitting Kit



		Buckle Closure	Hook & Loop
XS	ø 49 - 64mm	AS001B-PG	AS001-PG
S	ø 51 - 68mm	AS002B-PG	AS002-PG
M	ø 56 - 75 mm	AS003B-PG	AS003-PG
L	ø 75 - 94mm	AS004B-PG	AS004-PG
XL	ø 94 - 110mm	AS005B-PG	AS005-PG

Provide a soft supportive cuff that helps to secure the rear of the foot and ankle, while still allowing some movement of the foot.



All the belts and harnesses supplied with Plastic Buckles are **machine washable** at 60°c

FITTING KITS

All PostureGuard belts and harnesses are supplied with the Tri Glide system as standard, providing a simple, adjustable and versatile solution. To allow maximum flexibility, a full range of optional fitting kits is available to order separately, depending on individual requirements. This approach helps maintain a more competitive price point and reduces unnecessary waste, as additional fittings are not always required.

Tri Glide (included)

L9-P05-01A



Frame Clamp Fitting Kit

MC123-PG

Fits:

19mm

22mm

25mm



Locking Clamps



S

25mm WC25-PG

L

40mm WC40-PG

Triple Slot Mounting Plates



S

25mm MP012-PG

L

40mm MP034-PG

▶ AQUA RANGE

The Aqua Range is our **selection** of PostureGuard products also suitable for **wet environments**, supporting safe and effective postural management during hydrotherapy, bathing and personal care.

Introducing dedicated mounting kits for wet room applications, manufactured from corrosion resistant materials to ensure durability in wet environments.

	Padding (Individual)	Webbing & Buckle
S	150 x 40mm	25mm
M	170 x 60mm	40mm
L	230 x 70mm	40mm

All lower straps are 25mm wide

Pelvic Belt with Plastic Buckle

Padded | Dual Pull | Two Point

BPB-004-S-PG | BPB-004-M-PG | BPB-004-L-PG



Pelvic Belt with Plastic Buckle

Padded | Side Pull | Two Point

BPB-013-S-PG



Pelvic Belt with Plastic Buckle

Padded | Dual Pull | Four Point

BPB-005-S-PG | BPB-005-M-PG | BPB-005-L-PG



Pelvic Belt with Plastic Buckle

Padded | Side Pull | Four Point

BPB-014-S-PG



Support Strap | Buckle and Pull Tight Fastening



XS 28 x 8cm TS-001-PG



S 38 x 10cm TS-002-PG



M 47 x 13cm TS-003-PG




L 57 x 16cm TS-004-PG

Tri-Glide Fitting Kit included

Dynamic Harness

Quick Release, 25mm



S H009-PG **M** H010-PG **L** H011-PG

	A	B	C	D
S	110	180	160	45
M	145	230	200	50
L	175	260	230	60

All units are displayed in mm



Frame Clamp

MC123-AQ-PG

Fits:

19mm

22mm

25mm



Locking Clamps (pair)

WC25-AQ-PG | WC40-AQ-PG



Tri Glide (included)

L9-P05-01A



Fitting kits designed specifically for wet room environments, providing secure and corrosion resistant mounting for the products in the Aqua Range.

Why RMS?

- **With over 50 years of experience**, RMS has earned the trust of clinicians and the **NHS** through reliable, evidence-based products that improve lives.
- **Certified for quality and safety**, we operate under ISO 13485, with RMS products registered in accordance with the UK MDR and MHRA requirements, demonstrating compliance and our commitment to patient safety and product integrity.
- **End-user focus has been a core value since the beginning**, guiding every design to deliver comfort, confidence and independence through genuine understanding of individual needs.
- **Proudly built in Britain**, each product reflects our commitment to precision engineering, consistent quality and dependable local manufacturing.
- **A complete solution provider**, RMS supports clinicians from assessment to long-term care, delivering integrated posture and mobility solutions.
- **Nationwide coverage** with regional product specialists working alongside clinicians at assessments, handovers and reviews.



Rehabilitation Manufacturing Services

Thompson House, Unit 10, Styles Close,
Sittingbourne, Kent ME10 3BF

Tel: +44 (0) 1795 477280

Email: sales@rms-rehab.co.uk

www.rms-rehab.co.uk