



Buddy

CLASSROOM CHAIR



For over 50 years, RMS (Rehabilitation Manufacturing Services) has been designing and manufacturing postural support and mobility equipment for children and adults with complex needs.

Trusted by NHS Wheelchair Services, therapists and healthcare professionals in schools and CDCs across the UK, we have built a reputation for dependable products, clinical expertise and long-term support.



A TRUSTED PARTNER FOR SCHOOLS AND CDCS

As a British manufacturer, we maintain full control over design and production, ensuring consistent quality, reliable supply and the flexibility to meet individual requirements. Our range of seating systems, mobility products and postural support accessories is developed to improve comfort, participation and functional outcomes across educational, clinical and home environments.

Through our nationwide team of Product Specialists, we work closely with schools, therapists and families to assess needs, recommend appropriate solutions and provide ongoing support. From initial assessment through to delivery and follow-up, our aim is to make the process as straightforward and effective as possible for both staff and users.

WHY RMS

- **50 years** of specialist manufacturing experience
- **Trusted by NHS and Schools** across the UK
- **British-made**, quality-controlled production
- **Nationwide** assessment and product support
- **ISO 13485 certified** for quality and safety

Quality & Value

Without Compromise



Buddy

CLASSROOM CHAIR

Key Features

- **Durable** polypropylene construction
- **Highly adjustable**
- **Six sizes:** from early years to adulthood
- **Interchangeable** seat cushions
- Wide range of **accessories**
- **Value** for money
- **Highly recyclable**
- **Choice of colours**
- **Easy to clean** and maintain



Designed for Learning

Built to Last

The Buddy Classroom Chair is a durable and affordable seating solution designed for users requiring mild levels of postural support. Developed for use in schools, therapy environments and at home, the Buddy combines practical adjustability, a wide range of accessories and robust construction in a cost-effective package.



Manufactured from durable black polypropylene, the Buddy is built to withstand everyday use whilst remaining easy to clean. Quick-release hook-and-loop fastenings allow seat cushions to be removed easily for cleaning or exchanged with spare cushion sets between users, helping maintain a hygienic seating environment.

It is particularly well suited to busy educational and multi-user environments and is highly recyclable.

SIZED TO FIT THE CLASSROOM

The Buddy has been developed to integrate seamlessly into classrooms, therapy spaces and home environments.

Available in six sizes and compatible with furniture manufactured to BS EN 1729 standards, it provides a practical seating solution that helps users participate comfortably alongside their peers.

MOBILITY TO SUIT YOUR NEEDS

The Buddy is available with either a static base or a mobile base to suit different environments.

The mobile version features two fixed silent front wheels and two rear swivel castors with lockable brakes, providing smooth manoeuvrability and easy positioning.



STATIC BASE



WHEELED BASE

Size 1


Size 2

Size 3

Size 4

Size 5

Size 6



Seat Height	Seat Depth	Seat Width	Backrest Width	Backrest Height	Max User Weight
21.5 - 26.5 cm	19 - 24.5 cm	24.5 cm	22 cm	20.5 cm	25 kg
26 - 31.5 cm	25.5 - 31.5 cm	30 cm	27.5 cm	29.5 cm	45 kg
31 - 39 cm	30 - 36.5 cm	34.5 cm	32 cm	35 cm	50 kg
33.5 - 44 cm	34 - 41.5 cm	37.5 cm	35 cm	40 cm	60 kg
38.5 - 47.5 cm	40 - 48.5 cm	43 cm	40.5 cm	47 cm	70 kg
43 - 53 cm	46 - 55 cm	47.5 cm	45 cm	54 cm	80 kg

THOUGHTFUL COLOUR CHOICES

The Buddy is available in four carefully considered colours that balance modern classroom aesthetics with colours children naturally respond to. By moving away from a clinical appearance, the Buddy helps create a seating solution that children are happy to use as part of their everyday routine.



HOW TO ORDER

We recommend ordering the Buddy using our intuitive interactive prescription form to ensure the most appropriate configuration is selected.

For clinicians already familiar with the product, the Buddy can also be ordered directly through our website.



FLEXIBLE ACCESSORY OPTIONS

A wide range of accessories allows the Buddy to be adapted to suit individual requirements. Integrated fixing points throughout the chair enable the addition of accessories such as headrests, armrests, chest straps, wheels, brakes and push handles, providing additional support and functionality where required.



Flat Cushion



Ramped Cushion



Armrest Infill Pads



Black Tray



Flip-Up Footrest



Stability Skis



Push Handle



Flip Down Pommel



Type A



Headrests



Type D



Type G



Rollercoaster



Belts



Postural Support



Chest Straps



Pelvic Stabiliser



Harnesses



Manufactured in the UK

- **ISO 13485 certified**
- **Durable components** designed for long-term use
- Responsive **technical support** and aftersales assistance



The Buddy seating system is designed to be robust and provide safe, effective postural support for an average lifetime of 5-years, under normal use conditions.

Every component on the Buddy can be replaced to allow for continued use and re-issue following servicing & maintenance.

The actual service life of individual products may vary depending on frequency of use, environment and maintenance. The service life can be maximised by adhering to the manufacturer's IFU.

Disclaimer: *The above applies if the product is under normal use conditions, free from unreasonable or deliberate misuse.*



Scan to
download
the Buddy
Instructions

BOOK A FREE ASSESSMENT

Our nationwide team of Product Specialists can provide **assessments, demonstrations** and **product support**, helping schools configure each Buddy Seating System to meet the individual needs of the user.



Rehabilitation Manufacturing Services

Thompson House, Unit 10, Styles Close,
Sittingbourne, Kent ME10 3BF

Tel: +44 (0) 1795 477280

Email: sales@rms-rehab.co.uk

www.rms-rehab.co.uk